

Optogenetics applied to the High Throughput Screening: novel cellular assays controlled by light


Viviana Agus, Alberto Di Silvio, Jean-Francois Rolland, Anna Mondini,
Sara Tremolada, Katharina Montag, Lia Scarabottolo, Loredana Redaelli, Stefan Lohmer. AXXAM S.p.A, Bresso, Italy

Introduction

The recently developed Optogenetic technologies are based on genetically encoded elements which activity is controlled by light. By combining genetic targeting with optical stimulation, optogenetic tools can deliver a temporally precise stimulation or measurement of a biological process, in genetically defined group of cells or cell compartments.


Both "Light actuators" (proteins which transduce optical signals into physiological signals) and "Light sensors" (proteins which translate cell physiological signals into optical signals) have been expressed heterologously in different vertebrate systems, showing very good functionality even within intact tissues or behaving animals.

Due to the fact that they can be easily targetable to a specific cell type or cell compartment, and can be precisely and reversibly modulated by light pulses,



they have been largely used to study complex cellular mechanisms, such as electrical networking of neuronal cells, or in general to use the light stimulation to induce or measure a biological event.

We investigated the possibility of using these optogenetic tools to create a new generation of cellular assays suitable for the High Throughput Screening, where the light is used to modulate and measure the functionality of a gene of interest in a more precise, reversible and physiological way compared to the use of chemicals. In particular we verified if these recombinant optogenetic assays can be adapted to a miniaturized 384-well plate format, and also if the optics of a standard HTS instrumentation like the FLIPR (FLuorescence Imaging Plate Reader) could efficiently induce or detect the light modulation of these optogenetic tools.





Case study 1

Chr2 ion channel to modulate $Ca_v1.3$

In the first case study, we could demonstrate that the light-activated cation channel Chr2 D156A (from *Chlamydomonas reinhardtii*) expressed in HEK-293 cells can be very efficiently stimulated by the blue LEDs of the FLIPR instrument; the membrane depolarization signal obtained upon blue light stimulation was very strong and sustained, and could be very well detected in the miniaturized 384-well plate format used. By the co-expression of the $Ca_v1.3$ voltage-gated calcium channel, we could also demonstrate that the membrane depolarization induced by the channelrhodopsin stimulation is sufficient to elicit a very strong calcium influx through the $Ca_v1.3$ channel; we could also setup a protocol where, upon repetitive blue light stimulations, we

could study the $Ca_v1.3$ channel either in a "resting" or in an "inactivated" state. The optogenetic assay for $Ca_v1.3$ channel showed a very robust performance, and also pharmacological data in very good agreement with literature and with electrophysiology experiments.





Case study 2

bPAC enzyme to modulate HCN2

In the second case study, we stably expressed a light-activated enzyme, the Photoactivated Adenylyl Cyclase from *Beggiatoa* bacterium (bPAC) in HEK-293 cells.

By the co-expression with the HCN2 cyclic nucleotide gated channel, we could demonstrate that the blue LEDs of the FLIPR instrument can activate in a very fast and robust way the cAMP production through bPAC enzyme, resulting in a strong and sustained cation influx through the HCN2 channel.

The optogenetic assay for HCN2 channel showed an improved response window compared to the classical Forskolin stimulation of endogenous adenylyl cyclases, and pharmacological data in good agreement with literature.





Case study 3

ArcLight and RCaMP for simultaneous voltage and calcium imaging

In the third case study we co-expressed the Genetically Encoded Voltage Indicator ArcLight (from *Ciona intestinalis*) with the Genetically Encoded Calcium Indicator RCaMP in HEK-293 cells, and we demonstrated that this cell line can efficiently detect membrane voltage and intracellular calcium changes induced by

different ion channel targets. Both the green fluorescent ArcLight and red fluorescent RCaMP sensors showed a robust and homogeneous fluorescent signal, which can be very well detected by the FLIPR instrument LEDs, in a miniaturized 384-well plate assay format. Thanks to the nice spectral separation of these two sensors, we could

develop a protocol for the simultaneous detection of both membrane voltage and intracellular calcium changes in the same well, generating a multi-color, multi-parametric HTS-compatible cellular assay with improved content and reduced costs compared to the use of organic dyes.



Concluding remarks

- A new generation of recombinant cellular assays have been developed, by the use of light activated proteins and genetically encoded calcium or voltage indicators.
- The light-activated cation channel Channelrhodopsin 2 was proven to efficiently modulate the membrane voltage and the $Ca_v1.3$ calcium channel activity, allowing a robust and reliable detection of $Ca_v1.3$ state-dependent blockers.
- The light-activated adenylyl cyclase bPAC was proven to efficiently activate the HCN2 cyclic nucleotide gated cation channel, through the production of cAMP second messenger. The induced membrane depolarization was very fast and robust, showing higher performance compared to the classical Forskolin stimulation and good pharmacological behaviour.
- The Genetically Encoded Voltage Indicator ArcLight was proven to be able to detect the membrane voltage changes induced by a potassium, a sodium or a chloride ion channel target.
- The Genetically Encoded Calcium Indicator RCaMP was proven to be able to detect the intracellular calcium changes induced by a muscarinic receptor.
- The co-expression of ArcLight and RCaMP indicators can be used to generate a multi-color, multi-parametric reporter cell line suitable for the simultaneous detection of membrane voltage and intracellular calcium changes in the same well.
- Our results demonstrated that different optogenetic tools can be efficiently modulated by the FLIPR instrument optics, in a miniaturized 384-well plate format.
- The possibility to use optogenetic tools in an High Throughput Screening-compatible instrumentation, might deliver significant improvements in the development of future functional assays, allowing a temporal and spatial precision of stimulation and detection that would be superior to any other chemical or electrical methodology.

