

Selection of the most appropriate technology to develop an HTS grade SLC assay

Lia Scarabottolo, Sara Tremolada, Francesca Sassone, Loredana Redaelli

Axxam S.p.A., OpenZone - via A. Meucci, 3 - 20091 Bresso (Milan, Italy)

SLC assay development for screening purpose

The development of SLC functional assays for running High-throughput screening campaigns represents a fundamental approach to discover novel inhibitors and activators of SLC transporters for therapeutic purposes.

Different screening strategies have been used for the discovery of SLC modulators, including cellular substrate uptake assays with radioactivity or fluorescence detection; downstream functional assays with methods based on detection of pH and membrane potential; ligand binding assays with detection methods including radioactivity, fluorescence, mass spectrometry, or thermal shift (CETSA); phenotypic screening assays with different target identification methods.





In vitro, cell-free methods can also be used, anyway due to the SLC complex topography and numerous transmembrane domains, SLCs often require conditions that mimic their native environment to remain properly structured and functionally active; that's way cellular systems are preferable for the development of assays for screening purpose.

In the context of Resolute WP3 and WP6, Axxam has contributed, by mean of its consolidated experience and expertise, in the development of functional assays for many different SLCs, applying multiple and diversify technologies including fluorescent dyes, genetically encoded sensors, fluorescent substrate-based detection methods.

Assay requisites for screening purpose

To be suited for running a High Throughput Screening (HTS), an assay, for a cell line recombinantly expressing the SLC of interest, must satisfy some fundamental quality parameters:

- Possibility of miniaturization at least in 384 well plate format
- Specific response to the substrate (native or surrogate)
- Negligible response in mock transfected or non-induced cells
- Significant signal to background window, displaying a proper RZ' factor, which considers the signal amplitude and the standard deviation
- Signal stability over time
- Pharmacology reproducibility over time



Some examples

In this poster, some case studies of cell-based SLC assays, suitable for running HTS and developed by using different tools are shortly shown.

The detection systems used include Membrane Potential (MP) dye, SLC fluorescent substrate, Chloride and pH biosensors, Thallium dye, to record Thallium influx, as a surrogate for Potassium influx.

