

SLC12A2 assay development

Sara Tremolada, Viviana Agus, Loredana Redaelli, Lia Scarabottolo.

Axxam S.p.A., OpenZone - via A. Meucci, 3 - 20091 Bresso (Milan, Italy)

Introduction

- **SLC12A2** (coding for the electroneutral transporter **NKCC1**) mediates uptake of **Na⁺**, **K⁺**, and **Cl⁻** in a ratio of **1:1:2** and is inhibited by "loop diuretics" such as **furosemide** and **bumetanide**
- It is localized on the **plasma membrane (PM)**; it plays a role in **neuronal Cl⁻ homeostasis** and represents a **target for brain pathologies**
- It has also been shown to participate in the **maintenance and regulation of cell volume** (transport of ions is also associated to water transport)
- **Activation** of the transporter has been demonstrated by osmotic stress and is **regulated by phosphorylation**




Figure 1 - Structure of NKCC1 with important sites for Cl⁻ ions and Bumetanide binding, as well as phosphorylation sites.

- Treatment with **Bumetanide** has shown positive outcomes in clinical case studies of patients with **neurological/psychiatric conditions**. However, Bumetanide has a strong diuretic effect due to its inhibition of the kidney Cl⁻ transporter NKCC2 (SLC12A1) → critical issues with drug compliance and health concerns
- **Inhibitor compounds** that present **selectivity for NKCC1 versus NKCC2** may have a strong and consistent therapeutic value





Figure 2 - A novel NKCC1 inhibitor devoid of diuretic effects (ARN23746), may be clinically relevant for treatment of several neurological conditions characterized by impaired Cl⁻ homeostasis.


Assay development and HTS optimization was successfully carried out in 384-well plate format:

- **Cell line:** **HEK-293 Jump In T-REx/SLC12A2** (OE cell line from CeMM, N-terminally tagged SLC; IF expression confirmed, SLC localized on PM)
- **Read-out:** **Thallium (TI⁺) influx measurements** using a TI⁺-sensitive fluorescent dye (**Potassium assay kit**)

Thallium Dose-Response (D/R) and Bumetanide D/R




- Cells induction: **24h** (Doxycycline 1 µg/mL)
- Fluorescent Dye loading: **1h at RT** (in Chloride-free buffer)
- Double injection protocol (DMSO 0.5%, **FLIPRTETRA**). TI⁺ is injected in a buffer containing 25 mM Cl⁻




- TI⁺ D/R: significant difference between **induced and not-induced** condition
- Fluorescent signal is **fully abolished** by Bumetanide

DMSO sensitivity



DMSO up to 1.5% is not affecting:

- **EC₅₀/IC₅₀ values**
- **Assay window**





Multiplate test

A multiplate test was carried out on **5 plates** to assess the robustness of the assay and its eligibility for HTS purposes


- **MAX:** Bumetanide IC100 (30 µM, DMSO 0.5%)
- **MIN:** DMSO 0.5%
- **rZ' factor ≥ 0.5**
- Variability intraplate: ≤ **15%**
- Variability interplate: ≤ **20%**

Plate	RZ'
Plate 1	0.80
Plate 2	0.84
Plate 3	0.79
Plate 4	0.86
Plate 5	0.81



Assay quality criteria are fulfilled

Cell line stability



Cell line is stable up to p17:

- **EC₅₀ and IC₅₀ values are conserved**
- **Assay window is reproducible**

