

Live cell imaging of the Zinc solute carrier SLC30A8 (ZnT8)

Lucia Azzollini^{1#}, Dolores Del Prete^{1#}, Gernot Wolf², Christoph Klimek², Mattia Saggiaro¹, Fernanda Ricci¹, Eirini Christodoulaki², Alvaro Ingles-Prieto² & Lia Scarabottolo¹

¹ Axxam S.p.A., OpenZone, Bresso (Milan, Italy)

² CeMM Research Center for Molecular Medicine of the Austrian Academy of Sciences, Vienna, Austria

These authors equally contributed to this work

Introduction

- In pancreatic β -cells, zinc is crucial for the synthesis and secretion of insulin, which plays a key role in glucose homeostasis and which deficiency is the cause of diabetes.
- The accumulation of zinc in pancreatic cells is regulated by the solute carrier transporter SLC30A8 (Zinc Transporter 8, ZnT8), which transports zinc from cytoplasm in intracellular vesicles.
- Allelic variants of SLC30A8 gene have been linked to diabetes (in particular seem to increase the risk of the development of type 2 diabetes).





Figure 1 - Solute carrier family 30 (zinc transporter), member 8.

Challenges in developing an assay for SLC30A8	Solutions
SLC30A8 physiological expression in intracellular vesicles	SLC30A8 recombinant over-expression in plasma membrane
Cytotoxicity of SLC30A8 over-expression	Inducible system and optimization of induction condition
High background (endogenous Zinc transporters)	Imaging approach (single-cell analysis, high sensitivity)

A novel functional imaging- based assay specific for SLC30A8 was successfully developed in 384-well plate format:

- Cell lines: HEK293-JI cells overexpressing SLC30A8, SLC30A8 D110N_D224N, SLC30A8 R325W or R325Q (from CeMM, HA-tagged SLC; IF expression confirmed)
- Read-out: Zinc (Zn^{2+}) influx measurements using a Zinc-sensitive fluorescent dye (FluoZin-3, AM, cell permeant, by ThermoFisher)

Immunofluorescence for cellular localization of SLC30A8 over-expressed in HEK293 Jump-In T-REx cell lines




SLC30A8 D110N_D224N displayed the same expression levels and localization at the plasma membrane as SLC30A8.

The localization at the plasma membrane and overall expression of the diabetes variants (R325W or R325Q) was lower compared to SLC30A8 and the transport-dead mutant (D110N_D224N).

Figure 2 - Merge images of staining with: HA antibody for SLC30A8 and mutants (green), SLC1A5 antibody for plasma membrane (red), nuclei (blue). Images were generated by CeMM.

Image acquisition and image analysis



Cells stained with the FluoZin-3 (green) and counterstained with DAPI for nuclei (blue) were acquired by means of High Content Microscope (Opera Phenix, Revvity). Then cells were segmented to identify the nuclei using the DAPI channel.

Responding cells were selected based on a fluorescence intensity of FluoZin-3. The number of responding cells and relative mean fluorescence intensity were exported from the analysis.




Figure 3 - Image analysis process: 1. Nuclei segmentation with DAPI channel; 2. Cytoplasm ring region extended from nuclei segmentation; 3. Selection on responding cells; 4. Data extraction. Image analysis is performed with Harmony 4.9 (Revvity). Image acquired by 20x magnification.

Live imaging analysis of SLC30A8





Figure 4 - Representative images of FluoZin-3 staining after stimulation of induced cells with 1 mM $ZnCl_2$ (upper side); dose-response curves expressed as percentage of Zinc positive cells (lower side). Image acquired by 20x magnification.

- Cells induction: 24h (10 ng/mL doxycycline).
- Fluorescent dye loading: 45 mins at 37°C (in Ca^{2+} - Mg^{2+} -phosphate-free HBSS w 5mM EDTA).
- Automatic washing (3X) w Ca^{2+} - Mg^{2+} -phosphate-free HBSS.
- Stimulation w Zn^{2+} dose-response (as Zinc Chloride) up to 5 mM.
- Reading at Operetta CLS (Revvity).

• Zinc internalization was increased in SLC30A8 compared to unrelated transporter (SLC12A6) or WT cell line.

Comparison between $ZnCl_2$ transport activity across SLC30A8 mutants




Figure 5 - Representative images of FluoZin-3 staining after stimulation of induced cells with 1 mM $ZnCl_2$ (upper side); dose-response curves expressed as percentage of Zinc positive cells (lower side on the left); histograms representing percentage of Zinc positive cells at 1mM $ZnCl_2$. Statistical analysis: Brown-Forsythe and Welch ANOVA tests, p-values <0.05 = *, <0.005=**, <0.001=*** (lower side on the left). Image acquired by 20x magnification.

- Cytosolic accumulation of $ZnCl_2$ is greatly impaired in SLC30A8 D110N_D224N mutant and significantly lower than SLC30A8.
- Mutations associated with T2DM (SLC30A8 R325W and R325Q) display alterations in $ZnCl_2$ transport (higher n° of cells accumulating $ZnCl_2$ compared to SLC30A8).
- This cellular-based assay represents a novel and valuable tool to screen modulators of SLC30A8 transporter activity.